

SYNCHRON CX® Clinical Systems Operating Instructions Manuals

Important Notes for CX PRO Instruments Version 6.0

Important Notes

Purpose

This **Important Note** page consists of information or corrections that were discovered after the manual went to print. **These notes only apply to CX PRO instruments using software version 6.0.** In future manual revisions, this information will be incorporated into the body of the operations manual. After reading through the notes, place these **Important Notes** pages at the front of your operation manual for reference.

Note 1

Using Controls from Storage at -15° to -20° C:

When a Level # 3 control is removed from the freezer and immediately run on the instrument, the high viscosity of the control material may trigger a clot detection. Beckman Coulter Level #3 controls include *Decision*, *SYNCHRON*, *Triad* and *Nysspath*.

Note 2

Using a Floppy Diskette:

If you attempt to use a floppy diskette, a message (or similar message) "unable to mount the disk" may appear. The solution is to remove the floppy diskette, reboot the system, and use the diskette again.

Note 3

Defining QC Controls Sequentially:

To facilitate proper diagnostics and the transfer of information through the modem, QC control definitions must be created without gaps or blanks between each control item during QC setup.

Note 4

Performing a Data Backup/Restore Procedure:

When performing a backup/restore procedure (either a full backup or alignment data only), wait until the light on the floppy diskette drive turns off before removing the floppy diskette. This may take several seconds after the "Backup/Restore Complete" note appears on the screen. If the floppy diskette is removed before the light on the floppy drive turns off, the data transfer may be incomplete.